



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

ZymoPURE™ II Plasmid Maxiprep Kit

Catalog Nos. **D4202 & D4203** (Patent Pending)

Highlights

- Fast, easy, reliable, ultra-pure transfection grade plasmid DNA Maxiprep using a microcentrifuge spin-column.
- Innovative EndoZero™ technology enables transfection in sensitive cells and *in vivo* research (≤ 0.025 EU/ μ g of plasmid DNA).
- State-of-the-art ZymoPURE binding technology guarantees highly concentrated plasmid DNA directly from a spin-column in ≤ 20 minutes. No gravity filtration or ethanol precipitation!

Contents

| | |
|-----------------------------|------|
| Product Contents | 1 |
| Product Specifications..... | 1 |
| Product Description..... | 2 |
| Procedure Overview..... | 3 |
| Buffer Preparation | 4 |
| Protocols..... | 4-6 |
| Troubleshooting | 7 |
| Ordering Information | 8 |
| Related Products..... | 9-10 |

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

For Research Use Only

Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product, please call 1-888-882-9682.

Version 1.2.7

Product Contents:

| ZymoPURE™ II Plasmid Maxiprep Kit (Kit Size) | D4202 (10 preps.) | D4203 (20 preps.) | Storage Temperature |
|---|-------------------|-------------------|---------------------|
| ZymoPURE™ P1 ¹ (Red) | 150 ml | 2x 150 ml | 4°C |
| ZymoPURE™ P2 ² (Green) | 150 ml | 2x 150 ml | Room Temp. |
| ZymoPURE™ P3 (Yellow) | 150 ml | 2x 150 ml | Room Temp. |
| ZymoPURE™ Binding Buffer | 150 ml | 2x 150 ml | Room Temp. |
| ZymoPURE™ Wash 1 | 55 ml | 2x 55 ml | Room Temp. |
| ZymoPURE™ Wash 2 (Concentrate) | 23 ml | 2x 23 ml | Room Temp. |
| ZymoPURE™ Elution Buffer | 6 ml | 12 ml | Room Temp. |
| Zymo-Spin™ V-P Column Assemblies ³ | 10 | 20 | Room Temp. |
| ZymoPURE™ Syringe Filters | 10 | 20 | Room Temp. |
| ZymoPURE™ Syringe Plungers | 10 | 20 | Room Temp. |
| EndoZero™ Spin-Columns | 10 | 20 | Room Temp. |
| Collection Tubes | 10 | 20 | Room Temp. |
| Instruction Manual | 1 | 1 | - |

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

¹ ZymoPURE™ P1 contains RNase A (100 µg/ml) and is stable at room temperature without loss in RNase activity, however, for long-term storage the product should be stored at 4-8° C.

² Caution: ZymoPURE™ P2 Buffer contains NaOH. Please use proper safety precautions.

³ The Zymo-Spin™ V-P, 15 ml Conical Reservoir and 50 ml Reservoir are pre-assembled as a single unit.

Specifications:

- **DNA Purity:** Eluted DNA is ultrapure, endotoxin-free, and well suited for transfection, transformation, sequencing, restriction endonuclease digestion, *in vitro* transcription, *in vivo* studies, and other sensitive applications.
 - Typical Abs_{260/280} ≥ 1.8 and Abs_{260/230} ≥ 2.0
 - Endotoxin levels: ≤ 1 EU/µg of plasmid DNA using the Standard Protocol
≤ 0.025 EU/µg of plasmid DNA with optional EndoZero™ Spin-Column
- **Plasmid DNA Yield:** Up to 1.2 mg per preparation (*Actual yield is dependent on the plasmid copy number, culture growth conditions, and strain of E. coli utilized*)
- **Plasmid DNA Size:** Up to 200 kb
- **Recovery Volume:** ≥ 200 µl of ZymoPURE™ Elution Buffer or DNase-free water
- **Required Equipment:** Microcentrifuge and vacuum/vacuum manifold (recommended) or swinging bucket centrifuge.
- **Processing Time:** 20 min

Notes:

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

™ Trademarks of Zymo Research Corporation.

Several ZymoPURE™ product technologies are subject to U.S. and foreign patents or are patent pending.

pGL3™ is a registered trademark of Promega Corporation.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

Product Description


The **ZymoPURE™ II Plasmid Maxiprep Kit** features a simple spin-column based method for the purification of up to 1.2 mg of transfection grade plasmid DNA in less than 20 minutes. The eluted plasmid DNA is EndoZero and ready for immediate use in the most sensitive applications. The unique ZymoPURE methodology removes the need for slow gravity flow anion-exchange columns, alcohol precipitations, lengthy endotoxin removal incubations, and time-consuming centrifugation steps.

ZymoPURE™ technology uses a modified alkaline lysis method and features novel binding chemistry, which enables the highest yields and concentration of plasmid DNA (up to 3 µg/µl) directly from a spin-column. Coupling ZymoPURE with the innovative **EndoZero™ Spin-Columns**, to eliminate endotoxins, achieves EndoZero plasmid DNA (≤ 0.025 EU/µg of plasmid DNA), making it suitable for transfection, restriction endonuclease digestion, *in vivo* studies, bacterial transformation, PCR amplification, DNA sequencing, and other sensitive downstream applications.

As an added convenience, the **ZymoPURE™ II Plasmid Maxiprep Kit** contains colored buffers that permit error-free visualization and identification of complete bacterial cell lysis and neutralization. Syringe filters are included for rapid clearing of the lysate and the unique spin-column design allows the binding step to be performed using a vacuum or table top centrifuge.


For **Technical Assistance**, please contact **Zymo** at 1-888-882-9682 or E-mail tech@zymoresearch.com.

Highest Recovery




Plasmid DNA concentration and yield from the ZymoPURE II Maxiprep kit compared to two separate kits from Supplier Q. Plasmid DNA (pGL3[®]) was isolated from 150 ml of JM109 *E. coli* culture grown overnight following the manufacturer's suggested protocol (in duplicate). One (1) µl of eluted plasmid DNA was visualized post agarose gel electrophoresis. M, ZR 1 kb DNA Marker (Zymo Research).

Lowest Endotoxin Levels



Manufacturers' stated endotoxin for two separate Anion-Exchange kits from Supplier Q compared to ZymoPURE II.

Simplest Workflow



EZ-Load™
No Gravity Flow!




EZ-Elute™
No Alcohol
Precipitation!



EndoZero™
Endotoxins < FDA
limit for Vaccines

Procedure Overview:




Bacterial cells are resuspended in **ZymoPURE™ P1** (red).




The solution will turn dark purple and viscous following the addition of **ZymoPURE™ P2** (green) indicating bacterial lysis is complete.



The solution will turn yellow and a precipitate will form after adding **ZymoPURE™ P3** (yellow) indicating neutralization is complete.



The neutralized lysate is loaded into the **ZymoPURE™ Syringe Filter** and clarified into a new 50 ml conical tube.



ZymoPURE™ Binding Buffer is added to the cleared lysate and mixed thoroughly.



The mixture is loaded into the **Zymo-Spin™ V-P Column** using a vacuum manifold.



The **50 ml Reservoir** is removed and the **Zymo-Spin™ V-P Column** is washed using a vacuum manifold.



Ultra-pure plasmid DNA is eluted from the **Zymo-Spin™ V-P Column** using a microcentrifuge.



The eluted plasmid DNA is passed through the **EndoZero™ Column** using a microcentrifuge.

Buffer Preparation:

- ✓ Add 88 ml of 95% ethanol to the 23 ml **ZymoPURE™ Wash 2 (Concentrate)** before use.
- ✓ The **ZymoPURE™ P2** and **ZymoPURE™ Binding Buffer** may have precipitated. If this occurs, dissolve the precipitate by incubating the bottles at 30-37 °C for 10-20 minutes and mix by inversion. Do not microwave!

Before Starting:

- ✓ Centrifuge up to 150 ml of bacterial culture at $\geq 3,400 \times g$ for 10 minutes to pellet the cells¹. Discard supernatant.

Protocol:

The following procedure should be performed at room temperature (15-30°C).

1. Add 14 ml of **ZymoPURE™ P1 (Red)** to the bacterial cell pellet and resuspend completely by vortexing or pipetting.
2. Add 14 ml of **ZymoPURE™ P2 (Green)** and immediately mix by gently inverting the tube 6 times. Do not vortex! Let sit at room temperature for 2-3 minutes².

Cells are completely lysed when the solution appears clear, purple, and viscous.

3. Add 14 ml of **ZymoPURE™ P3 (Yellow)** and mix gently but thoroughly by inversion. Do not vortex!

The sample will turn yellow when the neutralization is complete and a yellowish precipitate will form.

4. Ensure the plug is attached to the Luer Lock at the bottom of the **ZymoPURE™ Syringe Filter**. Place the syringe filter upright in a tube rack and load the lysate into the ZymoPURE™ Syringe Filter³ and wait 5-8 minutes for the precipitate to float to the top.
5. Remove the Luer Lock plug from the bottom of the syringe and place it into a clean 50 ml conical tube. Place the plunger in the syringe and push the solution through the ZymoPURE™ Syringe Filter in one continuous motion until approximately 33-35 ml of cleared lysate is recovered. Save the cleared lysate!
6. Add 14 ml **ZymoPURE™ Binding Buffer** to the cleared lysate from step 5 and mix thoroughly by inverting the capped tube 10 times.

To continue processing the lysate using the recommended vacuum protocol, proceed to the next page. If a vacuum is not available, proceed to page 6 for an alternative centrifugation method.

Notes:

¹ A vessel with a minimum volume of 50 ml is required to prepare the bacterial lysate.

² Do not allow the lysis reaction to proceed for more than 3 minutes. Excessive lysis can result in denatured plasmid DNA.

³ If the precipitate has formed a homogenous layer at the surface of the neutralized lysate then invert the tube 3-4 times prior to loading the lysate into the **Zymo PURE™ Syringe Filter**.

Notes:

¹ To achieve optimal performance, the vacuum pump should be able to apply at least 400 mm Hg pressure. If less pressure is applied, centrifuge the column prior to washing to remove any residual lysate remaining in the matrix.

² The **ZymoPURE™ Elution Buffer** contains 10 mM Tris-HCl, pH 8.5 & 0.1 mM EDTA. If required, pure water can also be used to elute the DNA.

³ The DNA yield can be increased by pre-warming the **ZymoPURE™ Elution Buffer** to 50 °C and/or increasing the incubation period up to 10 minutes prior to centrifugation.

⁴ For low copy number plasmids or if higher concentration is desired, the plasmid DNA can be eluted in as little as 200 µl.

⁵ This optional step will reduce endotoxin levels from ≤ 1 EU/µg of plasmid DNA to ≤ 0.025 EU/µg of plasmid DNA.

Vacuum Protocol: (Recommended)

This product is compatible with any conventional vacuum-based manifold. The vacuum pump should be a single or double-staged unit capable of producing up to 400 mm Hg pressure at the vacuum manifold¹.

7. Ensure the connections of the **Zymo-Spin™ V-P Column Assembly** are finger-tight and place onto a vacuum manifold. (If vacuum is not available, see page 6 for the centrifugation protocol.)
8. With the vacuum off, add the entire mixture from step 6 into the Zymo-Spin™ V-P Column Assembly, and then turn on the vacuum¹ until all of the liquid has passed completely through the column.
9. Remove and discard the **50 ml Reservoir** from the top of the Zymo-Spin™ V-P Column Assembly.
10. With the vacuum off, add 5 ml of **ZymoPURE™ Wash 1** to the **15 ml Conical Reservoir**. Turn on the vacuum until all of the liquid has passed completely through the column.
11. With the vacuum off, add 5 ml of **ZymoPURE™ Wash 2** to the 15 ml Conical Reservoir. Turn on the vacuum until all of the liquid has passed completely through the column. Repeat this wash step.
12. Remove and discard the 15 ml Conical Reservoir and place the **Zymo-Spin™ V-P Column** in a **Collection Tube**. Centrifuge at $\geq 10,000 \times g$ for 1 minute, in a microcentrifuge, to remove any residual wash buffer.
13. Transfer the column into a clean 1.5 ml microcentrifuge tube and add 400 µl of **ZymoPURE™ Elution Buffer**^{2,3,4} directly to the column matrix. Wait 2 minutes, and then centrifuge at $\geq 10,000 \times g$ for 1 minute in a microcentrifuge.
14. *Optional:* For EndoZero Plasmid DNA⁵, remove the Luer Lock cap from the **EndoZero™ Spin-Column** and place the column in a clean 1.5 ml microcentrifuge tube. Add the entire eluate from step 13 into the EndoZero™ Spin-Column, wait 2 minutes, and then centrifuge at $5,000 \times g$ for 1 minute in a microcentrifuge. Store the eluted plasmid DNA at $\leq -20^{\circ}\text{C}$.

Centrifugation Protocol: (Alternative)

Perform steps 1-6 as indicated in the general protocol, see page 4.

7. Remove the **50 ml Reservoir** from the top of the **Zymo-Spin™ V-P Column Assembly**. Ensure the connection between the **15 ml Conical Reservoir** and **Zymo-Spin™ V-P column** is finger-tight and place the assembly into a 50 ml conical tube.
8. Add 14 ml of the mixture from step 6 into the **15 ml Conical Reservoir/Zymo-Spin™ V-P column assembly**, and then centrifuge the column at 500 x *g* for 2 minutes. Empty the 50 ml conical tube and repeat this step until the entire sample has passed through the column.
9. Add 5 ml of **ZymoPURE™ Wash 1** to the Zymo-Spin™ V-P column assembly and centrifuge the column at 500 x *g* for 2 minutes.
10. Add 5 ml of **ZymoPURE™ Wash 2** to the Zymo-Spin™ V-P column assembly and centrifuge the column for 2 minutes at 500 x *g*. Repeat the wash step.
11. Remove and discard the 15 ml Conical Reservoir and place the Zymo-Spin™ V-P Column into a **Collection Tube**. Centrifuge the column at $\geq 10,000$ x *g* for 1 minute, in a microcentrifuge, to remove any residual wash buffer.
12. Transfer the column into a clean 1.5 ml microcentrifuge tube and add 400 μ l of **ZymoPURE™ Elution Buffer**^{1,2,3} directly to the column matrix. Wait 2 minutes, and then centrifuge at $\geq 10,000$ x *g* for 1 minute in a microcentrifuge.
13. *Optional:* For EndoZero Plasmid DNA⁴, remove the Luer Lock cap from the **EndoZero™ Spin-Column** and place the column in a clean 1.5 ml tube. Add the entire eluate from step 12 into the EndoZero™ Spin-Column, wait 2 minutes, and then centrifuge at 5,000 x *g* for 1 minute in a microcentrifuge. Store the eluted plasmid DNA at $\leq -20^{\circ}\text{C}$.

Notes:

¹ The **Zymo PURE™ Elution Buffer** contains 10 mM Tris-HCl, pH 8.5, 0.1 mM EDTA. If required, pure water can also be used to elute the DNA.

² The DNA yield can be increased by pre-warming the **Zymo PURE™ Elution Buffer** to 50 °C and/or increasing the incubation period up to 10 minutes prior to centrifugation.

³ For low copy number plasmids or if higher concentration is desired, the plasmid DNA can be eluted in as little as 200 μ l.

⁴ This optional step will reduce endotoxin levels from ≤ 1 EU/ μ g of plasmid DNA to ≤ 0.025 EU/ μ g of plasmid DNA.

Troubleshooting Guide:

| Problem | Possible Causes and Suggested Solutions |
|---|--|
| Low DNA Yield | |
| <i>Culture growth conditions</i> | <ul style="list-style-type: none"> • Poor aeration of culture. The optimal culture volume to air volume ratio is 1:5 or less. For best aeration, use baffled culture flasks, or a vented or gas-permeable seal on the culture vessel. • The culture was overgrown, undergrown, contaminated, or antibiotics were omitted from the growth medium. Use a fresh culture for optimal performance. An OD₆₀₀ of 0.2-0.35 is the optimal optical density of a tenfold dilution of the culture. |
| <i>Cell density is too high</i> | <ul style="list-style-type: none"> • Too much culture used. Lysis and neutralization will be incomplete and the Zymo PURE™ Syringe Filter may clog during filtration. <u>More culture does not always equal more plasmid.</u> Incomplete lysis and neutralization are two of the most common causes of failed plasmid preps and both are caused by too much culture being used. • Incomplete lysis: After addition of ZymoPURE™ P2, the solution should change from opaque pink to a clear viscous purple, indicating complete lysis. Different <i>E. coli</i> strains often require different growth conditions and may vary in their susceptibility to alkaline lysis. • Incomplete neutralization: The solution should not be viscous following neutralization and the yellowish precipitate should appear fluffy and readily float to the surface. Make sure the neutralization is complete prior to filtration. Invert the tube an additional 2-3 times after the sample turns yellow following the addition of ZymoPURE™ P3. |
| <i>Lysate Clarification</i> | <ul style="list-style-type: none"> • Less than 33-35 ml of cleared lysate was recovered from the ZymoPURE™ Syringe Filter. For optimal performance, add 14 ml of ZymoPURE™ Binding Buffer to approximately 33-35 ml of clarified lysate. |
| <i>ZymoPURE P2 and ZymoPURE Binding Buffer precipitated</i> | <ul style="list-style-type: none"> • Both buffers may have precipitated during shipping. To completely resuspend the buffers, incubate the bottles at 30-37 °C for 10 minutes and mix by inversion. DO NOT MICROWAVE. |
| <i>Wash buffer</i> | <ul style="list-style-type: none"> • Ensure that ethanol has been added to the ZymoPURE™ Wash 2. • Ensure that the bottle cap is screwed on tightly after each use to prevent evaporation of the ethanol. |
| <i>DNA elution</i> | <ul style="list-style-type: none"> • Incomplete elution: For large size plasmids (> 10 kb), add ZymoPURE™ Elution Buffer and incubate the column for 5-10 minutes before centrifugation. Also, pre-warm the ZymoPURE™ Elution Buffer to 50 °C prior to elution. |
| Low DNA Quality | |
| <i>DNA does not perform well</i> | <ul style="list-style-type: none"> • Incomplete neutralization: Incomplete neutralization generates poor quality supernatant. Ensure that neutralization is complete by inverting the sample an additional 2-3 times after the addition of ZymoPURE™ P3 and extending the incubation. • Ethanol contamination in eluate. Centrifuge the Zymo-Spin™ V-P column as indicated in the protocol prior to adding the ZymoPURE™ Elution Buffer. |
| <i>RNA in eluate</i> | <ul style="list-style-type: none"> • Ensure that ZymoPURE™ P1 has been stored at 4°C. RNase A can be purchased separately if necessary. |
| <i>Genomic DNA in eluate</i> | <ul style="list-style-type: none"> • Improper handling (Sample was vortexed or handled too roughly). Genomic DNA contamination is usually caused by excessive mechanical shearing during the lysis and neutralization steps. Also, prolonged lysis or incomplete mixing of lysis or neutralization buffers may contribute to genomic DNA contamination in your sample. • Overgrown culture. Overgrown or old cultures may contain more genomic DNA contamination than fresh cultures. |

Ordering Information


| Product Description | Kit Size | Catalog No. |
|-----------------------------------|-----------|-------------|
| ZymoPURE™ II Plasmid Maxiprep Kit | 10 preps. | D4202 |
| ZymoPURE™ II Plasmid Maxiprep Kit | 20 preps. | D4203 |

| For Individual Sale | Amount | Catalog No. |
|---|--------|-------------|
| ZymoPURE™ P1 (Red) | 150 ml | D4200-1-150 |
| | 210 ml | D4200-1-210 |
| | 410 ml | D4200-1-410 |
| ZymoPURE™ P2 (Green) | 150 ml | D4200-2-150 |
| | 210 ml | D4200-2-210 |
| | 410 ml | D4200-2-410 |
| ZymoPURE™ P3 (Yellow) | 150 ml | D4200-3-150 |
| | 210 ml | D4200-3-210 |
| | 410 ml | D4200-3-410 |
| ZymoPURE™ Binding Buffer | 150 ml | D4200-4-150 |
| | 210 ml | D4200-4-210 |
| | 410 ml | D4200-4-410 |
| ZymoPURE™ Wash 1 | 55 ml | D4200-5-55 |
| | 410 ml | D4200-5-410 |
| ZymoPURE™ Wash 2 (Concentrate) | 23 ml | D4200-6-23 |
| ZymoPURE™ Elution Buffer | 6 ml | D4200-7-6 |
| | 12 ml | D4200-7-12 |
| | 30 ml | D4200-7-30 |
| Zymo-Spin™ V-P Column Assembly w/ 15 ml Conical and 50 ml Reservoir | 5 | C1042-5 |
| 15 ml Conical Reservoir | 25 | C1031-25 |
| 50 ml Reservoir | 25 | C1032-25 |
| ZymoPURE™ Syringe Filter | 5 | C1036-5 |
| ZymoPURE™ Syringe Plunger | 5 | C1037-5 |
| EndoZero™ Spin-Columns | 10 | C1051-10 |
| Collection Tubes | 50 | C1001-50 |
| | 500 | C1001-500 |
| | 1000 | C1001-1000 |

ZYMO RESEARCH CORP.

DNA PURIFICATION

What is Clean-Spin™ Technology?




The spin columns from Zymo Research have been designed to ensure complete elution with no binding/wash buffer carryover. The result is ultra-pure inhibitor-free DNA and RNA.

Purify DNA from PCR & other sources

DNA Clean & Concentrator™ (DCC™)

- ✓ Recovery of ultra-pure DNA that is free of salts and contaminants.
- ✓ Small ($\geq 6 \mu\text{l}$) elution volume.
- ✓ DNA is ideal for ligation, PCR, Next-Gen sequencing, etc.

| Product | Size (Cat. No.) |
|-----------------------------------|--|
| DNA Clean & Concentrator™-5 | 50 Preps. (D4013) 200 Preps. (D4014) |
| ZR-96 DNA Clean & Concentrator™-5 | 2 x 96 Preps. (D4023) 4 x 96 Preps. (D4024) |
| Genomic DNA Clean & Concentrator™ | 25 Preps. (D4010) 100 Preps. (D4011) |




High efficiency DNA recovery with the DCC™-5 compared to Supplier Q.

Boost DNA recoveries from agarose gels to >80%

Zymoclean™ Gel DNA Recovery

- ✓ Rapid (15 min.) recovery of ultra-pure DNA from agarose gels in $\geq 6 \mu\text{l}$.
- ✓ Ultra-pure DNA ideal for DNA ligation, sequencing, etc.
- ✓ Format also available for large DNA >20 kb.

| Product | Size (Cat. No.) |
|--|---|
| Zymoclean™ Gel DNA Recovery Kit | 50 Preps. (D4001) 200 Preps. (D4002) |
| Zymoclean™ Large Fragment DNA Recovery Kit | 25 Preps. (D4045) 100 Preps. (D4046) |



DNA fragments recovered from an agarose gel using the Zymoclean™ Gel DNA Recovery Kit. Lanes: M: DNA Ladder; 1-5: individual ladder DNA fragments.

Recover transfection-quality plasmid DNA directly from culture

Zyppy™ Plasmid Prep Kits

- ✓ The fastest, simplest method available for purifying high quality plasmid DNA from *E. coli*.
- ✓ Pellet-Free™ procedure omits conventional cell-pelleting and resuspension steps.
- ✓ Transfection quality plasmid DNA directly from culture in under 15 minutes.

Easy, Pellet-free Procedure: Add Lysis Buffer **Directly** to Bacterial Culture



| Product | Size (Cat. No.) |
|-----------------------------|--------------------|
| Zyppy™ Plasmid Miniprep Kit | 50 Preps. (D4036) |
| | 100 Preps. (D4019) |
| | 400 Preps. (D4020) |
| | 800 Preps. (D4037) |

OTHER INNOVATIVE PRODUCTS FROM ZYMO RESEARCH...


Competent cells for transformations *without* heat shock!

Mix & Go! Pre-made Competent *E. Coli*

- ✓ High efficiency: 10^8 - 10^9 transformants/ μ g plasmid DNA
- ✓ Just Mix & Go! Simply add DNA then spread. Transformation in as little as 20 seconds!

| Product | Size (Cat. No.) |
|--|---|
| Zymo 5 α (Same as DH5 α) | 10 x 100 μ l aliquots (T3007) |
| | 96 x 50 μ l aliquots (T3009) |
| | 96 x 50 μ l aliquots PCR-plate (T3010) |
| Zymo 10B (Same as DH10B) | 10 x 100 μ l aliquots (T3019) |
| | 96 x 50 μ l aliquots (T3020) |
| JM109 | 10 x 100 μ l aliquots (T3003) 96 x 50 μ l aliquots (T3005) |
| HB101 | 10 x 100 μ l aliquots (T3011) 96 x 50 μ l aliquots (T3013) |
| C600 | 10 x 100 μ l aliquots (T3015) |
| TG1 | 10 x 100 μ l aliquots (T3017) |

- ✓ No heat shock
- ✓ No incubations
- ✓ No outgrowth
- ✓ No wait!!!




* For Ampicillin selection only.


The fastest method for complete bisulfite conversion of DNA

EZ DNA Methylation-Lightning™ Kits

- ✓ The next generation of bisulfite conversion technology by the most cited provider in the industry
- ✓ Guarantees high conversion efficiencies of cytosine (>99.5%)
- ✓ Maintains the highest template integrity following bisulfite conversion
- ✓ Recovered DNA is ideal for PCR, MSP, array, bisulfite, and next-generation sequencing.



DNA Sequencing Results Following Bisulfite Treatment



Methylation Plot From Reduced Representation Bisulfite Sequencing (RRBS)

| Product | Size (Cat. No.) |
|--|-----------------------------------|
| EZ DNA Methylation-Lightning™ Kit | 50 rxns. (D5030) |
| | 200 rxns. (D5031) |
| EZ-96 DNA Methylation-Lightning™ Kit | Shallow-Well 2 x 96 rxns. (D5032) |
| | Deep-Well 2 x 96 rxns. (D5033) |
| EZ-96 DNA Methylation-Lightning™ MagPrep | 4 x 96 rxns. (D5046) |
| | 8 x 96 rxns. (D5047) |

ZYMO RESEARCH CORP.



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com