



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

# INSTRUCTION MANUAL

## **ZR-Duet<sup>TM</sup> DNA/RNA MiniPrep**

Catalog No. **D7001**

### **Highlights**

- Quick, 15 minute isolation and separation of DNA and RNA (~25 µg) from a wide range of sources using Clean-Spin<sup>TM</sup> column technology.
- High-quality DNA/RNA eluted in ≥25 µl is ready for reverse transcription, microarray, sequencing, etc.

### **Contents**

Product Contents .....	1
Specifications .....	1
Product Description .....	2
Buffer Preparation .....	3
Protocol .....	3, 4
Appendix .....	5
Ordering Information .....	6
Related Products .....	7

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please call 1-888-882-9682.

## Product Contents

<b>ZR-Duet™ DNA/RNA MiniPrep (Kit Size)</b>	<b>D7001 (50 preps.)</b>
<b>DNA/RNA Lysis Buffer</b>	50 ml
<b>DNA/RNA Prep Buffer</b>	50 ml
<b>DNA/RNA Wash Buffer<sup>1</sup> (concentrate)</b>	24 ml
<b>DNase/RNase-Free Water</b>	10 ml
<b>Zymo-Spin™ IIC Columns</b>	50
<b>Zymo-Spin™ IIIC Columns</b>	50
<b>Collection Tubes</b>	3x 50
<b>Instruction Manual</b>	1

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

**Storage Temperature** - Store all kit components (*i.e.*, buffers, columns) at room temperature.

<sup>1</sup> Ethanol must be added prior to use as indicated on **DNA/RNA Wash Buffer** label.

## Specifications

- **Sample Sources** – Cells, small amounts of *easy-to-lyse* tissue, buffy coat, buccal cells, plasma, serum, and other biological liquids. *Not compatible with whole blood.*<sup>2</sup>
- **Sample Size** – 10<sup>2</sup> to 5x10<sup>6</sup> cells in suspension or as tissue.
- **Recovery** – DNA and RNA can be eluted into small volumes (≥25 µl) allowing for a highly concentrated sample. Maximum DNA/RNA binding capacity of the provided columns is ~25 µg.
- **Size Limits** – Capable of recovering genomic DNA up to and above 40 kb. In most instances, mitochondrial DNA and viral DNA (if present) will also be recovered. Total RNA ≥17 nucleotides can be recovered.
- **Purity** – High quality genomic DNA and total RNA ( $A_{260}/A_{280} >1.8$ ,  $A_{260}/A_{230} >1.8$ ) is recovered. Traces of DNA may be present in the eluted RNA fraction. Trace DNA can be removed by DNase digestion (page 5).
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored at ≤-70 °C. The addition of RNase inhibitors is highly recommended for prolonged storage.
- **Equipment Needed** – Microcentrifuge

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility. RNAIate™ is a trademark of Ambion, Inc., Austin, Texas and is protected by various U.S. and foreign patents.


<sup>2</sup> For purification of DNA and RNA from whole blood, see the **ZR-Duet™ DNA/RNA MiniPrep Plus (D7003)**.

### ZYMO RESEARCH CORP.


## Product Description

The **ZR-Duet™ DNA/RNA MiniPrep** provides a quick method for the isolation of high quality genomic DNA and total RNA from small amounts of cells and tissue. The kit isolates *both* genomic DNA and a broad range of RNA species without the use of phenol. Small RNAs (*e.g.*, tRNAs, microRNAs) can be recovered following a simple adjustment within the RNA isolation protocol – *no extra steps are required!* Both DNA and RNA from up to  $5 \times 10^6$  cells can be eluted into volumes as little as 25  $\mu$ l in less than 15 minutes.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail [tech@zymoresearch.com](mailto:tech@zymoresearch.com).



Genomic DNA (lane 1, 2) and total RNA (lane 3, 4) isolated from human epithelial cells (HCT 116) with the **ZR-Duet™ DNA/RNA MiniPrep**.



PCR amplification of  $\beta$ -actin transcript (353 bp fragment shown) following DNA and RNA isolation from human epithelial cells (HCT 116) with the **ZR-Duet™ DNA/RNA MiniPrep**: PCR positive control (DNA template; lane 1, 2, 3), PCR negative control (RNA template; lane 4, 5, 6), RT-PCR (lane 7, 8, 9).

### ZYMO RESEARCH CORP.

Ensure the RNA isolation procedure is performed in an RNase-free environment.

**Notes:**

<sup>1</sup> In order to lyse samples completely, the amount of the **DNA/RNA Lysis Buffer** can be adjusted (*i.e.*, more buffer can be added).

<sup>2</sup> The capacity of the **Zymo-Spin™ IIC Column** is 800 µl. Columns can be reloaded to process volumes >800 µl.

<sup>3</sup> The maximum binding capacity of the **Zymo-Spin™ IIC and IIC Column** is ~25 µg of DNA/RNA.

## **Buffer Preparation**

Before starting, add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **DNA/RNA Wash Buffer** concentrate.

## **Protocol**

### **1. Sample Preparation**

- A. Adherent Cells:** Cells can be lysed directly in the culture container by removing liquid medium and adding **DNA/RNA Lysis Buffer**<sup>1</sup> directly to the monolayer (*e.g.*, 400 µl for 10<sup>2</sup> to 5x10<sup>6</sup> cells). Remove cells from the culture surface by pipetting, scraping, *etc.* Proceed to *Step 2*.
- B. Cells in Suspension:** Pellet the cells by gentle centrifugation (*e.g.*, 5 minutes at 500 x g). Remove the supernatant completely and resuspend the cell pellet in 400 µl **DNA/RNA Lysis Buffer**<sup>1</sup>. Vortex briefly. Proceed to *Step 2*.
- C. Solid Tissue Samples:** Add 400 µl **DNA/RNA Lysis Buffer**<sup>1</sup> to fresh or frozen tissue (up to ~25 mg) and homogenize the sample (*e.g.*, using a Dounce or similar homogenizer). Proceed to *Step 2*.
- D. Liquid Samples:** Add 3 volumes of **DNA/RNA Lysis Buffer**<sup>1</sup> for every volume of sample (*e.g.*, 300 µl buffer to 100 µl sample). Proceed to *Step 2*.

- 2. Transfer the sample from *Step 1* into a **Zymo-Spin™ IIC Column**<sup>2,3</sup> in the **Collection Tube** and centrifuge at ≥12,000 x g for 1 minute.**

**Save the flow-through for RNA purification and the column for DNA purification!**

**DNA Purification**

3. Transfer the **Zymo-Spin™ IIC Column** into a new **Collection Tube**.

4. Add 400 µl **DNA/RNA Prep Buffer** to the column and centrifuge at  $\geq 12,000 \times g$  for 1 minute. Discard the flow-through.
5. Add 700 µl **DNA/RNA Wash Buffer** to the column and centrifuge at  $\geq 12,000 \times g$  for 30 seconds. Discard the flow-through.
6. Add 400 µl **DNA/RNA Wash Buffer** and centrifuge the column for 2 minutes to ensure complete removal of the wash buffer. Carefully transfer the column into a clean microcentrifuge tube.

7. Add  $\geq 50$  µl **DNase/RNase-Free Water** directly to the column matrix and let stand 2 to 5 minutes at room temperature, then centrifuge at top speed for 30 seconds. The eluted DNA can be used immediately or stored at  $\leq -20^{\circ}\text{C}$ .

**RNA Purification**

3. Add 1 volume<sup>1</sup> ethanol (95-100%) to the flow-through in the **Collection Tube**<sup>2</sup> from *Step 2* and mix well by pipetting. Then transfer the sample into a **Zymo-Spin™ IIC Column** in a **Collection Tube** and centrifuge at  $\geq 12,000 \times g$  for 1 minute. Discard the flow-through.<sup>3</sup>

7. Add  $\geq 25$  µl **DNase/RNase-Free Water** directly to the column matrix and let stand 1 minute at room temperature. Centrifuge at  $10,000 \times g$  for 30 seconds. The eluted RNA can be used immediately or stored at  $\leq -70^{\circ}\text{C}$ .

**Notes:**

<sup>1</sup> Alternatively, to isolate RNAs  $\geq 200$  nt, add  $\frac{1}{2}$  volume ethanol (95-100%) to the sample flow-through.

<sup>2</sup> Capacity of the **Collection Tube** is 2 ml.

<sup>3</sup> At this point, RNA samples can be in-column DNase I treated (page 5).

**Notes:**

<sup>1</sup> Prior to use, reconstitute the lyophilized **DNase I** as indicated on the vial. Store frozen aliquots.

*\* Unit definition - one unit increases the absorbance of a high molecular weight DNA solution at a rate of 0.001 A<sub>260</sub> units/min/ml of reaction mixture at 25°C.*

**In-Column DNase I Digestion**

The DNase I digestion procedure can be performed using **DNase I Set** (E1010).<sup>1</sup>  
All centrifugation steps should be performed at 10,000 –16,000 x g for 30 seconds unless specified.

1. Wash the column with 400 µl **DNA/RNA Wash Buffer** and centrifuge. Discard the flow-through.

2. Add 80 µl **DNase I Reaction Mix** (below) directly to the column matrix.

<b>DNase I</b>	5 µl (1 U/µl)*
<b>DNA Digestion Buffer</b>	75 µl

3. Incubate the column at room temperature (20-30°C) for 15 minutes.  
Continue with RNA Purification: Page 4, Step 4.

**Ordering Information**

Product Description	Catalog No.	Kit Size
ZR-Duet™ DNA/RNA MiniPrep	D7001	50 Preps.

For Individual Sale	Catalog No.	Amount
DNA/RNA Lysis Buffer	D7001-1-50	50 ml
DNA/RNA Prep Buffer	D7010-2-10	10 ml
	D7010-2-25	25 ml
	D7010-2-50	50 ml
DNA/RNA Wash Buffer (concentrate)	D7010-3-6	6 ml
	D7010-3-12	12 ml
	D7010-3-24	24 ml
DNase/RNase-Free Water	W1001-1	1 ml
	W1001-4	4 ml
	W1001-6	6 ml
	W1001-10	10 ml
Zymo-Spin™ IIC Columns	C1011-50	50
	C1011-250	250
Zymo-Spin™ IIIC Columns	C1006-50	50
	C1006-250	250
Collection Tubes	C1001-50	50
	C1001-500	500
	C1001-1000	1000

**ZYMO RESEARCH CORP.**

## Related Products

Product	Description	Prep/Format	Catalog
<b>Total RNA Purification</b>			
ZR Whole-Blood RNA MiniPrep™	whole blood, partitioned blood	50/column 100/column	R1020 R1021
ZR-96 Whole-Blood RNA Kit™		2x96/plate	R1022
ZR Viral RNA Kit™	plasma, serum, liquids, cells, tissue	50/column 200/column	R1034 R1035
ZR-96 Viral RNA Kit™		2x96/plate 4x96/plate	R1040 R1041
ZR Urine RNA Isolation Kit™	urine, liquid samples	50/column	R1039
Quick-RNA™ MicroPrep	cells, tissue, buccal cells, buffy coat, plasma, serum, biological liquids	50/column	R1050
Quick-RNA™ MiniPrep		50/column 200/column	R1054 R1055
Quick-RNA™ MidiPrep		25/column	R1056
ZR-96 Quick-RNA™		2x96/plate 4x96/plate	R1052 R1053
ZR RNA MicroPrep™	cells, tissue, buccal cells, buffy coat, plasma, serum, biological liquids; DNA removal column, small-RNA recovery (≥17nt), <i>in-column</i> DNase treatment protocol	50/column 200/column	R1060 R1061
ZR RNA MiniPrep™		50/column 200/column	R1064 R1065
Pinpoint™ Slide RNA Isolation System Kit I	fresh tissue sections	50/column	R1003
Pinpoint™ Slide RNA Isolation System Kit II	paraffin-embedded tissue	50/column	R1007
ZR Fungal/Bacterial RNA MicroPrep™	bacteria, yeast, fungi; BashingBead™ lysis	50/column	R2010
ZR Fungal/Bacterial RNA MiniPrep™		50/column	R2014
ZR Plant RNA MiniPrep™	plant parts and tissues; BashingBead™ lysis, RT/PCR inhibitor removal	50/column	R2024
ZR Tissue & Insect RNA MicroPrep™	insect, small tissue samples; BashingBead™ lysis	50/column	R2030
YeaStar RNA Kit™	yeast, fungi	50/column	R1002
<b>RNA Clean-up, Concentration &amp; Recovery</b>			
RNA Clean & Concentrator™ -5	modified/labeled/impure/diluted RNA; small-RNA recovery (≥17nt); <i>acid phenol</i> extracted RNA directly from aqueous phase, <i>in-column</i> DNase treatment protocol	50/column 200/column	R1015 R1016
RNA Clean & Concentrator™ -25		50/column 100/column	R1017 R1018
RNA Clean & Concentrator™ -100		25/column	R1019
ZR-96 RNA Clean & Concentrator™		2x96/plate	R1080
DNA-Free RNA Kit™	DNase I treatment; small-RNA recovery (≥17nt)	50/column 200/column	R1013 R1014
Zymoclean™ Gel RNA Recovery Kit	agarose gel separated RNA	50/column	R1011
ZR small-RNA™ PAGE Recovery Kit	polyacrylamide gel separated RNA; small-RNA recovery (≥17nt)	20/column	R1070
<b>DNA/RNA Parallel Purification</b>			
ZR-Duet™ DNA/RNA MiniPrep	cells, tissue, liquids; DNA/RNA separation, small-RNA recovery (≥17nt), <i>in-column</i> DNase treatment protocol	50/column	D7001

### ZYMO RESEARCH CORP.