



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

ZR-96 Quick-RNA™

Catalog Nos. R1052 & R1053

Highlights

- High throughput (96-well) isolation of total RNA (including small RNAs) from a wide range of samples - single to 10⁶ cells.
- *DNA-free* RNA for use in any downstream application. *DNase I included.*

Contents

| | |
|----------------------------|------|
| Product Contents | 1 |
| Specifications | 1 |
| Product Description | 2 |
| Buffer Preparation | 3 |
| Protocols | 3, 4 |
| Appendices | 5 |
| Ordering Information | 6 |
| Related Products | 7 |

Satisfaction of all Zymo Research products is guaranteed. If you should be dissatisfied with this product please contact us.

Some difficult-to-lyse samples may require mechanical or enzymatic homogenization. For assistance, contact us at tech@zymoresearch.com.

Product Contents

| ZR-96 Quick-RNA™ (Kit Size) | R1052 (2x 96 Preps.) | R1053 (4x 96 Preps.) |
|--|---------------------------------|---------------------------------|
| RNA Lysis Buffer | 2x 100 ml | 4x 100 ml |
| RNA Prep Buffer | 100 ml | 2x 100 ml |
| RNA Wash Buffer¹ (concentrate) | 2x 48 ml | 4x 48 ml |
| DNase/RNase-Free Water | 10 ml | 30 ml |
| DNase I² (lyophilized) | 4 | 8 |
| DNA Digestion Buffer | 16 ml | 2x 16 ml |
| Silicon-A™ Plate | 2 | 4 |
| Collection Plate | 2 | 4 |
| Elution Plate | 2 | 4 |
| Cover Foil | 2 | 4 |
| Instruction Manual | 1 | 1 |

Note - Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Storage Temperature - Store all kit components (*i.e.*, buffers, columns) at room temperature.

¹ Before use, add 192 ml 100% ethanol (208 ml 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate.

² Prior to use, reconstitute the lyophilized **DNase I** as indicated on the vial prior to use. Store frozen aliquots.

Specifications

- **Sample Sources** – Cells or tissue samples, yeast, plant, bacteria, buccal cells, buffy coat, plasma, serum, and other biological liquids. *Compatible with DNA/RNA Shield™ and RNAlater™.*
- **Sample Storage** – Samples homogenized in RNA Lysis Buffer are stable and can be stored frozen prior to purification.
- **Sample Size** – Up to 10⁶ cells or 5 mg tissue.
- **RNA Purity** – High quality RNA ($A_{260}/A_{280} >1.8$, $A_{260}/A_{230} >1.8$) suitable for all downstream RNA-based manipulations.
- **RNA Recovery** – Up to 10 µg RNA can be eluted into ≥25 µl RNase-free water allowing for a highly concentrated sample.
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored frozen. RNase inhibitors can be included for prolonged storage.
- **Equipment Needed** – Centrifuge/rotor compatible with 96-well plates.

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility. RNAlater™ is a trademark of Ambion, Inc., Austin, Texas and is protected by various U.S. and foreign patents.


ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com


Product Description

The **ZR-96 Quick-RNA™** kit is an innovative product designed for the easy, reliable, and rapid isolation of DNA-free RNA from a wide range of cell (*up to 10⁶*) and tissue samples (*up to 5 mg*). The procedure combines a unique buffer system with Clean-Spin™ plate technology to yield high quality total RNA (*including small RNAs ~17-200 nt*) in about 30 minutes.

The procedure is simple: Add the provided **RNA Lysis Buffer** to a sample, then purify the RNA using the provided **Silicon-A™ Plate**. The result is highly-concentrated, *DNA-free* RNA that is suitable for subsequent RNA-based methods including RT-PCR, hybridization, sequencing *etc.*



The **ZR-96 Quick-RNA™** yields high quality total RNA. High levels of genomic DNA contamination are present in the preps from Suppliers Q & P but not with the **ZR-96 Quick-RNA™**. Total RNA was isolated from human epithelial cells (sans DNase treatment).




RNA isolated with the **ZR-96 Quick-RNA™** is DNA-free. Samples isolated with Supplier Q's kit are provided for comparison. Total RNA was isolated from 10⁶ human epithelial cells (with in-column DNase treatments for both kits). Each amplification curve represents an average of three independent isolation experiments.

For **Assistance**, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

Notes:

Use the **Direct-zol™-96 RNA** kit (Cat. Nos. R2054, R2055, R2056, R2057) for isolation of RNA *directly* (without phase separation) from samples in Trizol®, *etc.*

Use the **RNA Shield™** (Cat. Nos. R1100-50, R1100-250) for safe sample storage and transport at ambient temperatures.



The **Quick-RNA™** kits yield high quality RNA with high "RNA Integrity Numbers" (2200 TapeStation, Agilent).

ZYMO RESEARCH CORP.

Ensure the RNA isolation procedure is performed in an RNase-free environment.

Buffer Preparation

- ✓ Before starting, add 192 ml 100% ethanol (208 ml 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate.
- ✓ Reconstitute the lyophilized **DNase I** as indicated on the vial prior to use and store aliquots at -20°C.

Protocols

The RNA isolation consists of three steps: (I) *Sample Lysis/Homogenization*, (II) *Sample Clearing* and (III) *RNA Purification*.

All steps should be performed at room temperature (20-30 °C).

I. Sample Lysis/Homogenization

Notes:

Samples homogenized in **RNA Lysis Buffer** can be stored frozen for processing at a later time.

ZR Bashing Bead™ Lysis Tubes are available separately (Cat. Nos. S6002, S6003).

Processing plant tissue and other samples containing polyphenolics, humic acids, melanin, etc. may require use of the **OneStep™ PCR Inhibitor Removal Kit** (Cat. No. D6030).

Use the **DNA/RNA Shield™** for safe sample storage and transport at ambient temperatures.

Recommended **RNA Lysis Buffer** volumes

| RNA Lysis Buffer | 100 µl | 300 µl |
|-------------------------|-----------------------|-----------------------|
| Cells | Up to 10 ⁵ | Up to 10 ⁶ |
| Tissue | - | Up to 5 mg |

Adherent Cells

Lyse cells directly in the culture container by removing liquid medium and adding **RNA Lysis Buffer** directly to the monolayer.

Cells in Suspension

Pellet cells ($\leq 500 \times g$), remove the supernatant completely then resuspend the cell pellet in **RNA Lysis Buffer**. Vortex briefly.

Tissue and Tough-to-Lyse Samples

Fresh or frozen tissue (animal, plant, insect, yeast or bacteria) can be mechanically homogenized (e.g., **ZR BashingBead™ Lysis Tubes**) directly in the **RNA Lysis Buffer**.

Alternatively, tough-to-lyse tissue samples can be Proteinase K treated (page 5).

Liquids/Reaction Clean-up

DNase-treated RNA, labeling and *in vitro* transcription reactions can be processed directly by adding 4 volumes of **RNA Lysis Buffer** to each volume of sample (4:1) then mixing well.

Samples in DNA/RNA Shield™

Bring samples homogenized and stored in **DNA/RNA Shield™** to room temperature (20-30 °C). Then add 1 volume **RNA Lysis Buffer** (1:1), mix and proceed with Sample Clearing step.

Samples in DNA/RNA Shield™ can be Proteinase K treated (page 5).

Samples in RNA^{later}™

To process cells or liquids in **RNA^{later}™** (without reagent removal): Add 1 volume of RNase-free water or PBS to the sample (1:1). Then add 4 volumes **RNA Lysis Buffer** (4:1) and mix.

Alternatively, remove the **RNA^{later}™**, then proceed with Sample Lysis/Homogenization according to the sample type.

II. Sample Clearing

The following is recommended for cells and tissue (animal/plant) but can be omitted for cell-free liquids and low input samples ($\leq 10^5$ cells).

For particulate removal, centrifuge lysates at $\geq 12,000 \times g$ for 1 minute. Then transfer up to 300 μl of the supernatant into an RNase-free tube/plate (not provided).

III. RNA Purification

All centrifugation steps should be performed at $\geq 2,500 \times g$.

1. Add 1 volume ethanol (95-100%) to sample in the **RNA Lysis Buffer** [1:1] and mix well.
2. Transfer the mixture to a **Silicon-A™ Plate**¹ mounted on a **Collection Plate** and centrifuge for 5 minutes. Discard the flow-through.

3. In-column DNase I Treatment (optional)

This step can be used for trace DNA removal.

- a. Add 400 μl /well **RNA Wash Buffer** and centrifuge for 5 minutes. Discard the flow-through.
- b. For each sample to be treated, prepare **DNase I Reaction Mix** in an RNase-free tube (not provided). Mix well by gentle inversion:

| | |
|-----------------------------|------------------|
| DNase I ² | 5 μl |
| DNA Digestion Buffer | 35 μl |

- c. Add 40 μl **DNase I Reaction Mix** directly to the matrix. Incubate the plate at room temperature (20-30 °C) for 15 minutes. Then centrifuge for 5 minutes.

4. Add 400 μl **RNA Prep Buffer** to the plate and centrifuge for 5 minutes. Discard the flow-through.
5. Add 500 μl **RNA Wash Buffer** to the plate and centrifuge for 5 minutes. Discard the flow-through. **Repeat this step.**
6. Mount the **Silicon-A™ Plate** onto an **Elution Plate**, and add $\geq 25 \mu\text{l}$ **DNase/RNase-Free Water**³ directly to the matrix, then centrifuge for 5 minutes.

The eluted RNA can be used immediately or stored frozen. Use the **Cover Foil** to prevent evaporation.

Notes:

¹ To process samples $>600 \mu\text{l}$, **Silicon-A™ Plate** may be reloaded.

² Prior to use, reconstitute the lyophilized **DNase I** as indicated on the vial. Store frozen aliquots.

Unit definition - one unit increases the absorbance of a high molecular weight DNA solution at a rate of 0.001 A_{260} units/min/ml of reaction mixture at 25°C.

³ To maximize RNA yield, preheat the **DNase/RNase-Free Water** to 95° C, increase the elution volume and/or repeat the elution.

Proteinase K Digestion

Notes:

¹ **2X Digestion Buffer** (Cat. No. D3050-1-5 and D3050-1-20).

² **Proteinase K** (Cat. No. D3001-2-5 and D3001-2-20).

One unit of enzyme will hydrolyze urea-denatured hemoglobin to produce 1.0 μmole of tyrosine per minute at pH 7.5 at 37°C.

Example: up to 5 mg solid tissue or 10⁶ animal cells in DNA/RNA Shield™
2X Digestion Buffer¹
Proteinase K²

95 μl
95 μl
≥6 U

Prepare a Proteinase K reaction mix (see example above, scale-up as necessary). Incubate at 55°C for 30 minutes (*e.g.*, pelleted white blood cells) or 1-3 hours (solid tissue). Then add 1 volume **RNA Lysis Buffer** and proceed to Sample Clearing (page 4).

Ordering Information

| Product Description | Input | Binding | Catalog No. | Kit Size |
|-----------------------------|---|-------------|-------------|--------------|
| Quick-RNA™ MicroPrep | ~1-10 ⁶ cells | ~10 µg | R1050 | 50 Preps. |
| | | | R1051 | 200 Preps. |
| Quick-RNA™ MiniPrep | ~10 ² -10 ⁷ cells | ~100 µg | R1054 | 50 Preps. |
| | | | R1055 | 200 Preps. |
| Quick-RNA™ MidiPrep | ~10 ⁶ -10 ⁸ cells | ~1 mg | R1056 | 25 Preps. |
| ZR-96 Quick-RNA™ | ~1-10 ⁶ cells | ~10 µg/well | R1052 | 2x 96 Preps. |
| | | | R1053 | 4x 96 Preps. |

| For Individual Sale | Catalog No. | Amount |
|--------------------------------------|--|--------|
| RNA Lysis Buffer | R1060-1-50 | 50 ml |
| | R1060-1-100 | 100 ml |
| RNA Prep Buffer | R1060-2-10 | 10 ml |
| | R1060-2-25 | 25 ml |
| | R1060-2-100 | 100 ml |
| | R1003-3-6 | 6 ml |
| RNA Wash Buffer (concentrate) | R1003-3-12 | 12 ml |
| | R1003-3-24 | 24 ml |
| | R1003-3-48 | 48 ml |
| | DNase I (lyophilized) (250 U supplied with DNA Digestion Buffer, 4 ml) | E1010 |
| Silicon-A™ Plate | C2001 | 2 |
| Collection Plate | C2002 | 2 |
| Elution Plate | C2003 | 2 |
| Cover Foil | C2007-2 | 2 |
| | C2007-4 | 4 |
| DNase/RNase-Free Water | W1001-1 | 1 ml |
| | W1001-6 | 6 ml |
| | W1001-10 | 10 ml |

ZYMO RESEARCH CORP.

DNA PURIFICATION



What is Clean-Spin™ Technology?


The spin columns from Zymo Research have been designed to ensure complete elution with no binding/wash buffer carryover. The result is ultra-pure inhibitor-free DNA and RNA.

Purify DNA from PCR & other sources

DNA Clean & Concentrator™ (DCC™)

- ✓ Recovery of ultra-pure DNA that is free of salts and contaminants.
- ✓ Small ($\geq 6 \mu\text{l}$) elution volume.
- ✓ DNA is ideal for ligation, PCR, Next-Gen sequencing, etc.

| Product | Size (Cat. No.) |
|-----------------------------------|--|
| DNA Clean & Concentrator™-5 | 50 Preps. (D4013) 200 Preps. (D4014) |
| ZR-96 DNA Clean & Concentrator™-5 | 2 x 96 Preps. (D4023) 4 x 96 Preps. (D4024) |
| Genomic DNA Clean & Concentrator™ | 25 Preps. (D4010) 100 Preps. (D4011) |




High efficiency DNA recovery with the DCC™-5 compared to Supplier Q.

Boost DNA recoveries from agarose gels to >80%

Zymoclean™ Gel DNA Recovery

- ✓ Rapid (15 min.) recovery of ultra-pure DNA from agarose gels in $\geq 6 \mu\text{l}$.
- ✓ Ultra-pure DNA ideal for DNA ligation, sequencing, etc.
- ✓ Format also available for large DNA >20 kb.

| Product | Size (Cat. No.) |
|--|---|
| Zymoclean™ Gel DNA Recovery Kit | 50 Preps. (D4001) 200 Preps. (D4002) |
| Zymoclean™ Large Fragment DNA Recovery Kit | 25 Preps. (D4045) 100 Preps. (D4046) |



DNA fragments recovered from an agarose gel using the Zymoclean™ Gel DNA Recovery Kit. Lanes: M: DNA Ladder; 1-5: individual ladder DNA fragments.

Recover transfection-quality plasmid DNA directly from culture

Zyppy™ Plasmid Prep Kits

- ✓ The fastest, simplest method available for purifying high quality plasmid DNA from *E. coli*.
- ✓ Pellet-Free™ procedure omits conventional cell-pelleting and resuspension steps.
- ✓ Transfection quality plasmid DNA directly from culture in under 15 minutes.

Easy, Pellet-free Procedure: Add Lysis Buffer **Directly** to Bacterial Culture



| Product | Size (Cat. No.) |
|-----------------------------|--------------------|
| Zyppy™ Plasmid Miniprep Kit | 50 Preps. (D4036) |
| | 100 Preps. (D4019) |
| | 400 Preps. (D4020) |
| | 800 Preps. (D4037) |



RNA PURIFICATION

Get RNA *directly* from TRIzol® without phase separation


What is Clean-Spin™ Technology?
 The spin columns from Zymo Research have been designed to ensure complete elution with no binding/wash buffer carryover. The result is ultra-pure inhibitor-free DNA and RNA.

Direct-zol™ RNA

- ✓ For purification of high-quality small and large RNA *directly* from TRIzol®, TRI Reagent®, or similar.
- ✓ Bypasses phase separation and precipitation procedures allowing for unbiased recovery of miRNA

| Product | Size (Cat. No.) |
|--------------------------|---------------------|
| Direct-zol™ RNA MiniPrep | 50 Preps. (R2050) |
| | 50 Preps. (R2051)* |
| | 200 Preps. (R2052) |
| | 200 Preps. (R2053)* |

96-well and MagBead formats also available!
 DNase I included in all kits.
 * Supplied with TRI-Reagent®




High-quality small and large RNA are effectively recovered with the Direct-zol™ kit. RNA is DNA-free.

Isolate DNA-free RNA from 1 to 10⁷ cells in minutes

Quick-RNA™

- ✓ Isolation of total, large, or small RNA – *You decide!*
- ✓ Ultra clean, high-quality RNA from a single cell to 10⁷ cells.
- ✓ DNA-free RNA ideal for any downstream application – *DNase I included.*

| Product | Size (Cat. No.) |
|----------------------|-----------------------|
| Quick-RNA™ MicroPrep | 50 Preps. (R1050) |
| | 200 Preps. (R1051) |
| Quick-RNA™ MiniPrep | 50 Preps. (R1054) |
| | 200 Preps. (R1055) |
| ZR-96 Quick-RNA™ | 2 x 96 Preps. (R1052) |
| | 4 x 96 Preps. (R1053) |



Isolate total, large, or small RNA with the Quick-RNA™ kit. RNA is DNA-free using the Quick-RNA™ kit

Purify RNA from enzymatic and labeling reactions in 5 minutes

RNA Clean & Concentrator™

- ✓ Recover ultra-pure RNA in small (≥6 µl) elution volumes.
- ✓ Compatible with TRIzol®, phenol, chloroform, and RNase inhibitors (RNAlater®).
- ✓ RNA is ideal for RT-PCR, q-PCR, hybridization, arrays, RNA interference, etc.

| Product | Size (Cat. No.) |
|---------------------------------|--------------------------|
| RNA Clean & Concentrator™-5 | 50 Preps. (R1015) |
| | 200 Preps. (R1016) |
| RNA Clean & Concentrator™-25 | 50 Preps. (R1017) |
| | 100 Preps. (R1018) |
| ZR-96 RNA Clean & Concentrator™ | 2x96 well plates (R1080) |
| DNA-Free RNA Kit™ | 50 Preps. (R1013) |
| | 200 Preps. (R1014) |



The following are trademarks of other companies: pGEM®, Promega Corp.; TRIzol® and TRI Reagent®, Molecular Research Center, Inc.; DH5® and DH10B™, Life Technologies, Inc.



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com